



Cystoscopy Accessories

Cystoscope Sheath



Cystoscopy sheath are used to examine the inside of the bladder & urethra and urinary tract with use of Cystoscope.

Cystoscopy sheath comes with obturator.

| Features: | Sizes: |
|---|---------|
| Endoscopic imaging support. | 15.5 Fr |
| Tapered ends and assist in penetration. | 17 Fr |
| Durable, Safe and sterile instruments for each patient. | 19 Fr |
| Compatible with all cystoscopes . | 20 Fr |
| With 2 fixed stopcocks. | 22 Fr |
| Compatible with Karl STORZ. | 23.5 Fr |
| | 25 Fr |

Single Channel Cystoscope Bridge



Single Channel Cystoscope Bridge provides extra Instrument channel for Cystoscope Procedure.

Features:

Designed to provide excellent utility.

Durable & highly compatible.

Advanced configuration .

Autoclavable. .

Compatible with Karl STORZ and Richard Wolf

Double Channel Cystoscope Bridge



Double Channel Cystoscope Bridge provides two extra Instrument channel for Cystoscope Procedure.

Features:

Designed to provide excellent utility.

Durable & highly compatible.

Advanced configuration .

Autoclavable. .

Compatible with Karl STORZ and Richard Wolf

Cystoscopy Accessories

Optical Biopsy Forcep



Use for tissue identification by optical analysis and biopsy sampling at a site within the body.

Features:

Durable & highly reliable.


Improved quality parameters.

Anti-corrosive.

Fine functionality.

Compatible with Karl STORZ.

Flexible Biopsy Forcep



Use for tissue identification by optical analysis and biopsy sampling at a site within the body.

Features:

Superior cutting performance.

Less bleeding.

Autoclavable & Reusable.

High longevity.


Sizes

Length

5 Fr, 7 Fr

450 mm

Stent Removing Forceps



Stent Removing Forceps is designed for the easy removal of stents.

Features:

Reusable & Autoclavable.

Strong grasping.

Alligator jaw for strong grip.

Sizes

Length

5 Fr, 7 Fr

450 mm

Bladder Stone Removal Instruments

Stone Crushing Forcep



Stone Crushing Forcep is use to remove bladder stone. Compatible with 23.5 Fr or 25 Fr Cystoscopy Sheath.

Features:

Available : Stone Crushing Forcep

Provide excellent utility.

Advanced configuration.

Stone Crushing Forcep SET

Reusable & Autoclavable.

(Contains Stone Crushing Forcep and 25 Fr

Cystoscopy Sheath)

Compatible with Karl STORZ and Richard Wolf

Stone Punch SET / Stone Punch Urology



Stone Punch Urology SET use for remove bladder stones.

Features:

Fine finish & Corrosion resistance.

Longer life & Easy to use.

High performance.

Advanced configuration.


Stone Punch SET Includes:

Stone Punch

Sheath

Obturator

Litho Bridge



Litho Bridge is use with Pneumatic Lithotripter for bladder stone Fragmentation.

Features:

Reusable & Autoclavable.


Compatible with 23.5 Fr or 25 Fr Cystoscope Sheath.

High performance.


Compatible with Karl STORZ and Richard Wolf

Ureteroscopy Accessories


URS Forcep



URS Alligator Forcep



URS Tripronge Forcep



URS Flat Jaw Forcep

URS Forcep used for the retrieval of fragmented stones from ureteral and renal pelvis after performing Lithoclast procedure.

Features:

Durable & Easy to use.

Reusable & Autoclavable.

Anti-corrosive.

Smooth finish.

High performance.

Strong grasping.

Sizes:

3 Fr, 4 Fr, 5 Fr

Length:

610 mm

Types:


URS Alligator Forcep

URS Tripronge Forcep

URS Flat Jaw Forcep

Nephroscopy Accessories

PCNL Forcep



PCNL Alligator Forcep

PCNL Alligator Forcep Used for the retrieval of fragmented stones from kidney after performing Lithoclast procedure.

Features:

Reusable & Autoclavable.

Durable & Easy to use.

Anti-corrosive.

Smooth finish.

High performance.

Strong grasping.

Sizes:

6 Fr, 9 Fr, 12 Fr

Length:

360 mm


Types:

PCNL Alligator Forcep


PCNL Tripronge Forcep

PCNL Peanut Jaw Forcep


PCNL Flat Jaw U Handle Forcep



PCNL Tripronge Forcep



PCNL Peanut Jaw Forcep



PCNL Flat Jaw U Handle Forcep

Endoscopic Accessories

Alkine Telescopic Metal Dilator SET / PCN Metal Dilator SET



Alkine Telescopic Metal Dilator SET is used in dilation during PCNL Procedure.

Features:

Reusable & Autoclavable.

Optimum finish.

Reliable nature.

Anti-corrosive.

Sizes:

6 Fr - 30 Fr

Optical Urethrotome Sheath

Optical Urethrotome Sheath is used in Urethrotomy Procedure which involves incision of the urethra for relief of a stricture.

21 FR Optical Urethrotome Sheath comes with Obturator.

Features:

Good Locking System.

Tapered ends and assist in penetration.

With 2 fixed stopcocks.

Precise design.

Anti-corrosive.

Available with:

Optical Urethrotome Sheath

Optical Urethrotome Sheath with Halfmoon



Optical Urethrotome Sheath



Optical Urethrotome Sheath with Halfmoon

Endoscopic Accessories

OTIS Urethrotome



An Otis urethrotomy is the blind transurethral incision of the urethra by an Otis urethrotome.

Features:


Reusable & Autoclavable.

Anti-corrosive.

Smooth finish & Strong construction.

High performance & Precise design.

Bugbee Electrodes



Monopolar flexible bugbee electrodes are used in cystoscopes, resectoscopes or ureteroscopes for coagulation/ ablation of tissue.

Sizes:

3 Fr to 9 Fr

Hey groove dilator SET



Hey groove dilator set contains 4 no. dilators.

Sizes:

4/7, 8/11, 11/14, 14/16

Endoscopic Accessories

Male Urethral Dilator SET



Male urethral dilator set contains 12 Dilators

Sizes:

28/32, 26/30, 24/28, 22/26, 20/24, 18/22, 16/20, 14/18, 12/16, 10/14, 8/12, 6/10

Available:

Male urethral dilator set - With Gide Wire Passing

Male urethral dilator set - Without Gide Wire Passing

Female Urethral Meatal Dilator SET / Hegar Dilators




Hegar dilator SET contain 10 Dilators.

Sizes:

3/3.5, 4/4.5, 5/5.5, 6/6.5, 7/7.5, 8/8.5, 9/9.5, 10/10.5, 11/11.5, 12/12.5

Female Urethral Teflon Dilator SET



Female urethral Teflon dilator SET contain 5 Teflon Dilators.

Sizes:

13/14, 11/12, 9/10, 7/8, 5/6

Address:

Advin Urology (Division Of Advin Health Care)

424-Maruti Plaza, Sardar Patel Chowk Road, krishnanagar, Saijpur Bogha, Ahmedabad-382345, Gujarat, India.

E-Mail: exports@advinurology.com

Phone: +91-76007 27250

Web: www.advinurology.com